

EasiVac[®]

Disposable Prostate Chip Evacuator

Evacuate, analyze, dispose.

Easy Connection

- Adaptors to allow direct connection to various resectoscope sheaths
- Connection fills container and bulb with irrigation fluid to create suction
- Leakproof connection to prevent spills and reduce clean up

Suction and Filtration

- Rapidly collects fragments remaining after resection
- Anti-reflux system prevents chips from reentering into bladder
- Strong bulb allows for no wait between pressures for quick evacuation

Smart, Durable, and Disposable Construction

- Shock proof construction will not break in case of drop
- Main chamber can be closed with lid and sent to lab for analyzation
- Disposable and not made with natural rubber latex

For ordering and additional information please contact your sales representative or Coloplast Interventional Urology Surgical Support at **800-258-3476**



Ordering Information

Coloplast Interventional Urology Surgical Support **800-258-3476**

Ref Number	Description	UOM	Sterilization
E80103	Easivac® Bladder Evacuator with Wolf Connector	Box of 20	Sterile
E80203	Easivac® Bladder Evacuator with STORZ Connector	Box of 20	Sterile
E80303	Easivac® Bladder Evacuator with OLYMPUS Connector	Box of 20	Sterile
E80403	Easivac® Bladder Evacuator with CIRCON-ACMI Connector	Box of 20	Sterile

Indications for Use: Removal of prostate chips after TURP.

Warnings: Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

Adverse Events: The following events have been reported although their occurrence greatly depends on medical conditions of patient: infection, encrustation, obstruction, rupture, migration, bladder irritation symptoms, pain, hematuria, erosion.

See Instructions for Use for detailed information regarding warnings/precautions, adverse events prior to using this product. For further information contact Coloplast Corp at 1-800-258-3476 and/or consult the company website at www.coloplast.us.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

[Ostomy Care / Contenance Care / Wound & Skin Care / Interventional Urology](#)

Coloplast Corp. Minneapolis, MN 55411 / Interventional Urology Surgical Support 1-800-258-3476

www.coloplast.us The Coloplast logo is a registered trademark of Coloplast A/S.

© 2020 Coloplast Corp. All rights reserved.

