


Ureteral Catheters

Efficient urinary tract access, drainage, and navigation




Flush Ureteral Catheter

- Various tip designs to facilitate access to the urinary tract
- Suitable for use with a guidewire during endoscopic procedures
- Constructed of PEBA material

Product Description	Body Diameter	Length	Side Eyes	Guidewire Compatibility	Product Code
Flush Ureteral Catheter Straight Open Tip	4.8 Fr	74 cm	—	0.038"	ACP205
	6 Fr	74 cm	—	0.038"	ACP206
	7 Fr	74 cm	—	0.038"	ACP207
	8 Fr	74 cm	—	0.038"	ACP208
Flush Ureteral Catheter Open Tip with Side Eyes	3 Fr	74 cm	2	0.021"	ACP303
	4 Fr	74 cm	2	0.035"	ACP304
	4.8 Fr	74 cm	2	0.038"	ACP305
	6 Fr	74 cm	2	0.038"	ACP306
	7 Fr	74 cm	2	0.038"	ACP307
Flush Ureteral Catheter Coude Open Tip	3 Fr	74 cm	—	0.021"	ACP603
	4 Fr	74 cm	—	0.035"	ACP604
	4.8 Fr	74 cm	—	0.038"	ACP605
	6 Fr	74 cm	—	0.038"	ACP606
	7 Fr	74 cm	—	0.038"	ACP607
Flush Ureteral Catheter Coude Open Tip with Side Eyes	3 Fr	74 cm	2	0.021"	ACP403
	4 Fr	74 cm	2	0.035"	ACP404
	4.8 Fr	74 cm	2	0.038"	ACP405
	6 Fr	74 cm	2	0.038"	ACP406
	7 Fr	74 cm	2	0.038"	ACP407

Hydrophilic-coated Floppy Tip Ureteral Catheter


- Hydrophilic coating reduces coefficient of friction of the catheter
- Floppy tip for smooth passage through ureter
- Radiopaque tip for improved visualization



Product Description	Body Diameter	Length	Side Eyes	Guidewire Compatibility	Product Code
Hydrophilic-coated Ureteral Catheter Neoplex Floppy tip	5 Fr	70 cm	—	0.035"	ACR205
	6 Fr	70 cm	—	0.035"	ACR206
	7 Fr	70 cm	—	0.035"	ACR207

Chevassu Catheters for Retrograde UreteroPyelography (RUP)

- Allows for effective injection of contract medium
- Various body and cone/bulb diameter offerings
- Constructed of Neoplex material



Product Description	Body Diameter	Cone/Bulb Diameter	Length	Side Eyes	Guidewire Compatibility	Product Code
Chevassu Ureteral Catheter Compatible with Guidewire	5 Fr	7 Fr	73 cm	—	0.035"	AC5C07
Chevassu Ureteral Catheter	3 Fr	4 Fr	70 cm	—	—	AC5904
	3 Fr	5 Fr	70 cm	—	—	AC5905
	4 Fr	6 Fr	70 cm	—	—	AC5906
	5 Fr	7 Fr	70 cm	—	—	AC5907
	5 Fr	8 Fr	70 cm	—	—	AC5908
	6 Fr	9 Fr	70 cm	—	—	AC5909
	7 Fr	10 Fr	70 cm	—	—	AC5910
	7 Fr	11 Fr	70 cm	—	—	AC5911
	7 Fr	12 Fr	70 cm	—	—	AC5912

Ordering Information

For ordering and additional information please contact your sales representative or Coloplast Interventional Urology Surgical Support at **800-258-3476**

Indications for Use:

Interventional catheters (flush, open straight, open elbowed and open floppy tip catheters): Endourological procedures, to restore or maintain drainage of the upper urinary tract and to inject contrast medium or saline.
Catheters for Retrograde UreteroPyelography (Chevassu catheters): injection of contrast medium or saline during endourological procedures.

Warnings: Injecting contrast medium may cause side effects. It is therefore advisable to respect the use precautions specified in the instructions for the product used.

Adverse Events: Infection, encrustation, obstruction, migration, bladder irritation symptoms, and pain.

Ostomy Care / Continence Care / Wound & Skin Care / Interventional Urology

