

X-flow[®]

High performance silicone prostatic catheter

Performance and Safety

- Large bevel and eyes to ensure excellent drainage capacity.
- Wide irrigation and drainage channels to prevent blocking.
- Closed round tip for ease of placement of an introducer when needed.
- Stiff walls to facilitate clot evacuation.
- Clear catheter body to monitor drainage quality and visualize urine.
- All Silicone, to avoid any risk of allergy or toxicity.

A range of tips designed to meet patient needs

Dufour

Delinotte

Couvelaire

Straight

Phthalate and latex free

Ordering Information

Coloplast Urology Care Service Support **800-258-3476**

X-flow silicone prostatic catheters

Balloon Volume	Ways	Tip	18Fr	20Fr	22Fr	24Fr
50 ml	3	Dufour	AB6318	AB6320	AB6322	AB6324
		Couvelaire	AB6118	AB6120	AB6122	AB6124
		Straight	AB6018	AB6020	AB6022	AB6024
		Delinotte	AB6218	AB6220	AB6222	AB6224
	2	Dufour	AB6418	AB6420	AB6422	AB6424
		Couvelaire	AB6518	AB6520	AB6522	AB6524
		Straight	AB6618	AB6620	AB6622	AB6624
		Delinotte	AB6718	AB6720	AB6722	AB6724
30 ml	3	Dufour	AB6A18	AB6A20	AB6A22	AB6A24
		Straight	AB6C18	AB6C20	AB6C22	AB6C24

Indications for Use: For short-term drainage of bladder urine, postoperative bladder irrigation-lavage, and for haemostasis of the prostatic fossa following prostate surgery. .

Warnings:

- Using an iodine-based irrigation-lavage solution may damage some silicone catheters.
- Do not inflate the balloon beyond the maximum value indicated.
- Do not clamp the catheter. Use a plug if necessary.
- If the catheter needs to be secured, the adhesive must be applied to the connector.
- To lubricate catheters without a coating, a water-based lubricant is recommended.
- Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

Adverse Events: Several adverse events have been described with the use of balloon catheters. Some are related to the patient's conditions, the others to the procedure or the device: bladder irritation symptoms, pain, urinary tract infection, encrustation and stone formation, urinary tract trauma.

See Instructions for Use for detailed information regarding warnings/precautions, adverse events prior to using this product. For further information, contact Coloplast Corp at 1-800-258-3476 and/or consult the company website at www.coloplast.com.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

[Ostomy Care](#)
[Continence Care](#)
[Wound & Skin Care](#)
[Urology Care](#)

Coloplast Corp. Urology Care Surgical Support 800.258.3476
www.us.coloplast.com The Coloplast logo is a registered trademark of Coloplast A/S. © 2017. All rights reserved.

PM-01684 8.17

